

OPERATOR MANUAL

**Lab 250 and Lab 110**  
**Life Sciences Small Sterilizers**

(11/13/06)

P387352-364

**Rev. 0**




**Table 1. Recommended Loads by Sterilizer Size**

						Volume (Liters)
16	16	26"	(406	406	660 )	109
20	20	38"	(508	508	965 )	249

**Table 2. Sterilizer Configurations**

16	16	26"	(406	406	660 )	
20	20	38"	(508	508	965 )	
16	16	26"	(406	406	660 )	
20	20	38"	(508	508	965 )	


**Advisory**






5	CONTROL INTERFACE .....	5-1
5.1	.....	5-1
5.1.1	.....	5-1
5.1.2	.....	5-2
5.2	.....	5-4
5.2.1	.....	5-4
5.2.2	..... ( * )	5-7
5.2.3	..... ( * )	5-10
5.3	.....	5-13
5.4	..... ( * )	

7	ALARMS .....	7-1
8	ROUTINE MAINTENANCE .....	8-1
.1	.1.1	-1
.1	.1.2 /	-7
.1	.1.3	-
.1	.1.4	-
.2	.3	-11
.3	.4	-13
.4	.5	-14
.5	.6	-15
.6	.7	-15
9	PREVENTIVE MAINTENANCE .....	9-1
.1	.2	-1




WARNING-ELECTRIC SHOCK AND BURN HAZARD:



WARNING-PERSONAL INJURY HAZARD:



WARNING-BURN HAZARD:



---

WARNING-BURN HAZARD:









---

CAUTION - POSSIBLE EQUIPMENT DAMAGE:

- ⚠ Do not use the device in a wet or damp environment.
- ⚠ Do not use the device near open flames or heat sources.
- ⚠ Do not use the device in a confined space.

## Definition of Symbols

Symbol	Definition
	High voltage
	Ground
	Warning
	Warning
A	Amperes
V	Volts
~	AC
Hz	Hertz
$\phi$	Phase
SN	Signal-to-noise ratio



## 2.2 Technical Specifications

### 2.2.1 Electrical Service

120V AC, 60Hz, 15A, 1.0 phase.  
 (1.3 phase / 1.0 phase).  
 3

### 2.2.2 Sterilizer Final Check

1.0 phase / 1.0 phase.  
 1.0 phase / 1.0 phase.  
 1.0 phase / 1.0 phase.

### 2.2.3 Cycle Operation

1.0 phase / 1.0 phase.  
 (1.3 phase / 1.0 phase).  
 (121° 250°F).

### 2.2.4 Overall Size and Weight





<b>16 x 16 x 26"</b>	660	•	1	2	02	•	•
<b>Sterilizer:</b>	(26"	•	74.5"	•	35.5"	•	•)
<b>20 x 20 x 38"</b>	762	•	1	2	1152	•	•
<b>Sterilizer:</b>	(30"	•	74.5"	•	45.375"	•	•)

• **Electric:**


	230	1.5	1	50/60
	120	.5	1	50/60

2.2.7 Utility Requirements  
For Units Equipped With  
Optional Electric Steam  
Sterilizers


• **Electric:**

	120	2.5	1	50/60
	20	3.2	3	50/60
	440/480	31/37	3	50/60
	240	72.2	3	50/60

• **Hot Water:**


	137.140°F	344.7	(20, 50)
---	-----------	-------	----------

• **Stainless-Steel Option Only:**

	1	Ω	
---	---	---	--

2.2.8 Environmental  
Conditions

• **Cold Water:**

	206.70°F (21°C)	344.7	(30, 50)
---	-----------------	-------	----------

10° 32° (50° 0°F)

	10	0%	2	(5)
---	----	----	---	-----

# TECHNIQUES OF STERILIZATION

---

---

## 3.1 Recommended

Glassware 0025 Tc0 Tw9.864 Manual

3.-09(6ommended)Tj3.34.011 inverted,(3.2.16ommended)

**Table 3-2. Gravity Cycle Parameters (cont.)**

Minimum

Minimum


---

## 3.2 Control Measures for Verifying Sterilization Process

---

### 3.2.1 Biological Monitors

### 3.2.2 Testing for Prevacuum Efficiency

 **WARNING – STERILITY ASSURANCE HAZARD:** Load sterility may be compromised if the biological indicator or vacuum leak test indicates a potential problem. If these indicators show a potential problem, refer the situation to a qualified service technician before using the sterilizer further.



---

### 3.3 Bowie-Dick Test

---



---

### 3.4 Vacuum Leak Test

---

---

### 3.5 Recommendations for the Sterilization Process

---



Table 3-4

Table 3-4

### 3.7 Recommendations for Sterilizing Liquids

**⚠ WARNING – EXPLOSION HAZARD:** This sterilizer is not designed to process flammable compounds.

**⚠ WARNING – PERSONAL INJURY HAZARD:** Avoid personal injury from bursting bottles. Liquid sterilization cycle must only be used for liquids in borosilicate flasks with vented closures.

**⚠ WARNING – BURN HAZARD:** When sterilizing liquids, *always* observe the following procedures:

- Use Liquid cycle only.
- Use only vented closures.
- Use only Type I borosilicate glass bottles.
- Do not allow hot bottles to be jolted.

**⚠ WARNING – BURN HAZARD:** Steam may be released from the chamber when the door is opened. Step back from the sterilizer each time the door is opened to minimize contact with steam vapor.

**⚠ CAUTION:** Sterilization of chloride-containing solutions (e.g., saline) can cause chamber corrosion as discussed by the manufacturer. If, however, chloride-containing solutions must be processed, clean the chamber after each use.

IMPORTANT:




Figure 3-1. Vented Closures







```

===== GRAVITY =====
=====
CYCLE START AT      XX:XX:XX A/P
                   ON      XX/XX/XX

CYCLE COUNT          0
LOGIN NAME:          XXXXXX
STERILIZER           XXXXXX

CYCLE TYPE           GRAVITY
CYCLE NO.            2

STER TEMP = 121.0 C
CONTROL TEMP = 122.5 C
STER TIME = 0:30:00
DRY TIME = 0:01:00

- TIME              T=C          V=inHg
                   P=psig
-----
C 11:48:24A        66.7           0.3V
C 11:49:24A       112.7          10.0P
S 11:49:43A       121.2          16.6P
S 11:51:43A       122.6          17.7P
S 11:53:43A       123.3          17.8P
S 11:55:43A       123.6          16.8P
S 11:57:43A       122.6          17.0P
S 11:59:43A       122.6          17.2P
S 12:01:43P       122.5          17.0P
S 12:03:43P       122.4          17.2P
S 12:05:43P       122.5          16.8P
S 12:07:43P       122.4          16.9P
S 12:11:43P       122.4          17.0P
S 12:13:43P       122.5          17.1P
S 12:15:43P       122.6          17.0P
S 12:17:43P       122.7          17.0P
S 12:19:43P       122.6          16.8P
E 12:19:44P       122.6          16.9P
E 12:19:54P       113.7           3.2P
E 12:20:03P       99.9           11.1V
E 12:21:03P       40.5           28.1V
Z 12:21:46P       68.4            0.5V

LOAD                020903

CHAMBER TEMP MAX=124.8 C
CHAMBER TEMP MIN=121.2 C

CONDITON           = 1:19
STERILIZE          = 30:01
EXHAUST            = 1:42
TOTAL CYCLE        = 33:02

=====
===== READY TO UNLOAD =====
=====

```


Figure 4-5. Co de sed Pri tout

*(Faint musical notation or decorative elements)*

NOTE: E *(faint text)* a a/a a...




#### 4.2.3 Operating Mode



**Figure 4-6. Pri tout: Sterilizer Type a d Ma ufacturer**

#### 4.2.4 Cycle Start



NOTE: C o /a a u a , n S P r / o r M o .


#### 4.2.5 End-of-Cycle Performance Summary



#### 4.2.6 Alarm Condition



NOTE: R r o SECTION 7, ALARMS, o r , o p o u a a r o g o .



**Figure 4-7. Pri tout: Type of Alarm, Chamber Temperature a d Chamber Pressure, a d Time of Occurre ce**

---

## 4.3 Manual Operation of Manual Door

---

1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 23, 24, 25, 26, 27, 28, 29, 30, 31, 32, 33, 34, 35, 36, 37, 38, 39, 40, 41, 42, 43, 44, 45, 46, 47, 48, 49, 50, 51, 52, 53, 54, 55, 56, 57, 58, 59, 60, 61, 62, 63, 64, 65, 66, 67, 68, 69, 70, 71, 72, 73, 74, 75, 76, 77, 78, 79, 80, 81, 82, 83, 84, 85, 86, 87, 88, 89, 90, 91, 92, 93, 94, 95, 96, 97, 98, 99, 100





---

## 4.5 Optional Electric Steam Generator

---


٤٥

**Table 4-1. Required Feed Water Quality for Carbo Steel Steam Generators**


Condition	Normal Recommended	Maximum Recommended
		140°F (60°C)
3*	0-17 µS/cm	130 µS/cm
	50-150 µS/cm	250 µS/cm
3	50-100 µS/cm	100 µS/cm
	6.0 - 7.5	6.5 - 7.5
	0.1 - 1.0 µS/cm	2.5 µS/cm
	2000-6000	26,000



\*17.1 µS/cm - 100 µS/cm

1. CONTROL INTERFACE



Screen #84

2.  ( F# 4- ).



NOTE:  ,   
**CANCEL**

10. **CONTINUE**


The image shows a musical score for a piece titled "CONTINUE". The score is written on a single staff with a treble clef. The key signature has one sharp (F#), and the time signature is 4/4. The piece begins with a 10-measure introduction, indicated by the number "10." on the left. The word "CONTINUE" is written in large, bold, capital letters across the first few measures. The notation includes various rhythmic values such as quarter notes, eighth notes, and sixteenth notes, along with rests and dynamic markings. The piece concludes with a double bar line.



## 5.1 General Description



### 5.1.1 Information and Help



STERIS



Screen #22

EXIT


(; #22). EXIT

### 5.1.2 Password Entry


STERIS



Screen #11



Screen #11




INCORRECT PASSWORD. 11 X

## 5.2 Operating Modes


### 5.2.1 Gravity Sterilizer Only

OPERATING MODE



Screen #2

RIGHT ARROW




Screen #19

LEFT ARROW


DOOR OPEN, READY, JACKET CHARGE

CYCLE



Screen #17.5

RIGHT ARROW



Screen #17.5

LEFT ARROW

CANCEL CYCLE  
START CYCLE

NOTE: If a door is closed, the cycle will not start. The door must be closed before the cycle can start.


# PRINT CYCLE VALUES

Item	Quantity	Unit Price	Total Price
1.	1	2 00	2 00
		0 30 00	0 30 00
		121.0	121.0
		6.0	6.0
		1.5	1.5
		1.0	1.0
	2	10.0	20.0
		0 15 00	0 15 00
		0	0
-----			
2.	2	2 00	4 00
		0 30 00	0 30 00
		121.0	242.0
		6.0	12.0
		1.5	3.0
		1.0	2.0
	2	10.0	20.0
		0 15 00	0 15 00
		0	0
-----			
3.	3	2 00	6 00
		0 30 00	0 30 00
		121.0	363.0
		6.0	18.0
		1.5	4.5
		1.0	3.0
	2	10.0	20.0
		0 15 00	0 15 00
		0	0
-----			
4.	4	2 00	8 00
		0 30 00	0 30 00
		121.0	484.0
		6.0	24.0
		1.5	6.0
		1.0	4.0
	2	10.0	20.0

		10.0	20.0
		0 15 00	0 15 00
		0	0
-----			
5.	5	2 00	10 00
		0 30 00	0 30 00
		121.0	605.0
		6.0	30.0
		1.5	7.5
		1.0	5.0
	2	10.0	20.0
		0 15 00	0 15 00
		0	0
-----			
6.	6	2 00	12 00
		0 30 00	0 30 00
		121.0	726.0
		6.0	36.0
		1.5	9.0
		1.0	6.0
	2	10.0	20.0
		0 15 00	0 15 00
		0	0
-----			
7.	7	2 00	14 00
		0 45 00	0 45 00
		121.0	847.0
		6.0	42.0
		1.5	10.5
		1.0	7.0
	2	10.0	20.0
		0 15 00	0 15 00
		0	0
-----			
8.	8	2 00	16 00
		0 45 00	0 45 00
		121.0	968.0
		6.0	48.0
		1.5	12.0
		1.0	8.0
	2	10.0	20.0
		0 15 00	0 15 00
		0	0
-----			




STANDARD MODE  
 DOOR OPEN, READY, 115.0 (23.0 F), JACKET CHARGE  
 CYCLE



Screen #17.5

RIGHT ARROW



Screen #17.5

LEFT ARROW CANCEL CYCLE  
 START CYCLE

NOTE: CLOSE DOOR(S)




PRINT CYCLE VALUES

/ /	
1.	1
	2 00
	4
	10.0
	15.0
	0 30 00
	121.0
	6.0
	1.5
	1.0
	2
	10.0
	0 15 00
	0
-----	
2.	2
	2 00
	4
	10.0
	15.0
	0 30 00
	121.0
	6.0
	1.5
	1.0
	2
	10.0
	0 15 00
	0
-----	
3.	3
	2 00
	4
	10.0
	15.0
	0 30 00
	121.0
	6.0
	1.5
	1.0
	2
	10.0
	0 15 00
	0
-----	
4.	4
	2 00
	4
	10.0
	15.0
	0 30 00
	121.0
	6.0
	1.5
	1.0
	2

		10.0
		0 15 00
		0
-----		
5.	5	
		2 00
		0 30 00
		121.0
		6.0
		1.5
		1.0
		2
		10.0
		0 15 00
		0
-----		
6.	6	
		2 00
		0 30 00
		121.0
		6.0
		1.5
		1.0
		2
		10.0
		0 15 00
		0
-----		
7.	7	
		2 00
		0 45 00
		121.0
		6.0
		1.5
		1.0
		2
		10.0
		0
-----		
8.	8	
		2 00
		0 45 00
		121.0
		6.0
		1.5
		1.0
		2
		10.0
		0

# OPERATING MODE




**ISO MODE**  
 NOTE: I # a P a a a a # a 80.0 C (176.0 F), a # a  
 o r /



**STANDARD MODE**  
 DOOR OPEN, READY, 115.0 (23 .0 F).

RIGHT ARROW



Screen #17.5

LEFT ARROW

CANCEL CYCLE  
START CYCLE

NOTE: If a door is closed, the cycle will not start.

# PRINT CYCLE VALUES

			/ /
1.	1		
		2 00	
		0 30 00	
		121.0	
		6.0	
		1.5	
		1.0	
		2	
		10.0	
		0 15 00	
		0	
-----			
2.	2		
		2 00	
		0 30 00	
		121.0	
		6.0	
		1.5	
		1.0	
		2	
		10.0	
		0 15 00	
		0	
-----			
3.	3		
		2 00	
		0 45 00	
		121.0	
		6.0	
		1.5	
		1.0	
		2	
		10.0	
		0	
-----			
4.	4		
		2 00	
		0 45 00	

		1.0	
		2	
		0	
-----			
5.78			
		0 30 00	
		78.0	
		6.0	
		1.5	
		1.0	
		2	
-----			
6.88			
		0 30 00	
		88.0	
		6.0	
		1.5	
		1.0	
		2	
-----			
7.104			
		0 30 00	
		104.0	
		6.0	
		1.5	
		1.0	
		2	
-----			
6.104			
		0 30 00	
		104.0	
		6.0	
		1.5	
		1.0	
		2	
-----			

## 5.3 Status Buttons


PRES/TEMP    PRES/TEMP, MACHINE I/O, GRAPH  
 PRES/TEMP

NOTE: LOAD 1 a, F, a, n, o, a, p, p, o, r, t.



Screen #47

EXIT  
 MACHINE I/O



Screen #56

EXIT EXIT



Scree #2

PRINT STATUS .F

11 25 52 50.0 0.0

		م	م
		10.0	
		0 15 00	
		0	م
		-----	
1.	1	2 00	
		0 30 00	
		121.0	
		6.0	
		1.5	
		1.0	م
		2	م
		10.0	م
		0 15 00	
		0	م
		-----	
2.	2	2 00	
		0 30 00	
		121.0	
		6.0	
		1.5	
		1.0	م
		2	م
		10.0	م
		0 15 00	
		0	م
		-----	
3.	3	2 00	
		0 30 00	
		121.0	
		6.0	
		1.5	
		1.0	م
		2	م
		10.0	م
		0 15 00	
		0	م
		-----	
4.	4	2 00	
		0 30 00	
		121.0	
		6.0	
		1.5	
		1.0	م
		2	م

		10.0	
		0 15 00	
		0	م
		-----	
5.	5	2 00	
		0 30 00	
		121.0	
		6.0	
		1.5	
		1.0	م
		2	م
		0	م
		-----	
6.	6	2 00	
		0 30 00	
		121.0	
		6.0	
		1.5	
		1.0	م
		2	م
		0	م
		-----	
7.	7	2 00	
		0 45 00	
		121.0	
		6.0	
		1.5	
		1.0	م
		2	م
		0	م
		-----	
8.	8	2 00	
		0 45 00	
		121.0	
		6.0	
		1.5	
		1.0	م
		2	م
		0	م
		-----	
		م	م
		-----	
		1م	
		5م	
		60م	
		1م	
		30م	
		10م	
		60م	
4	4 0	1.	

5م





Scree #4





Scree #4







## 5.4.2 Printer Operation (Optional Mylox Printer)



---

## 5.5 Compact Flash Card

---


300,000 300,000 300,000 300,000

**VERY IMPORTANT:**

6 ( )




Screen #6



---

## 5.6 RS232 Download To PC



---






## 6.1 Before Operating Sterilizer


### CAUTION - POSSIBLE EQUIPMENT DAMAGE:


- Gasket must be fully retracted prior to opening sterilizer door.
- If 0 dry time is selected, sterilizer automatically initiates a vapor removal phase in place of dry cycle. This phase can still draw a vacuum to 5 in Hg. Consult device manufacturer's recommendations to ensure devices being processed can withstand this depth of vacuum.
- Lifting the chamber float switch when cleaning the chamber may cause the sterilizer control to initiate a Chamber Flooded alarm. If this alarm condition occurs, the operator must turn the control power OFF then ON to clear the alarm. The control power switch is located in the mechanical area at the side of the sterilizer. Placement of the sterilizer indicator does not clear this alarm.
- Allow thermostatic traps to cool down to room temperature before removing cover. Since there is nothing to limit expansion, the bellows may rupture or fatigue if trap is opened while hot.
- Actuation at less than 75% of rated pressure can allow debris to contact the seat and cause the safety valve to leak. A leaking safety valve must be replaced.
- Inefficient service clearance will make repairs more difficult and time-consuming.
- Piping sized too small may cause water hammer, resulting in damage to the sterilizer.
- After installation, it is mandatory to brace piping at the drain funnel so that it will not move vertically.
- Make sure door opening is clear of any obstruction before closing the door(s).
- Do not attempt to open sterilizer door during manual operation unless chamber is at 0 psig (0 bar).
- Never use a wire brush, abrasives, or steel wool on door and chamber assembly. Do not use cleaners containing chloride on stainless-steel surfaces. Chloride-based cleaners will deteriorate stainless steel, eventually leading to failure of the vessel.
- Immediately wipe up saline solution spills on loading car, to prevent damage to stainless steel.
- Do not use cleaners containing chlorides on loading cars. Chloride-based cleaners will deteriorate the loading car metal.
- Sterilization of chloride-containing solutions (e.g., saline) can cause chamber corrosion and is not recommended by the manufacturer. If, however, chloride-containing solutions must be processed, clean the chamber after each use.
- Avoid damage to the integral steam generator daily. Flush the generator daily. Failure to flush generator daily will void the manufacturer's warranty.

1.  ( F  6-1).

2.  M

3.  ( F  6-2).

4. 

5.  M



6.  $115^{\circ}$  (23 °F).


7. 6.1.1,

1. ( F 6-3).

2. M

3.

4.



**Figure 6-4: Chamber Door Half Way Out of Chamber**

1. ( F 6-4).





---

## 6.2 Gravity Cycle

---




### WARNING – BURN HAZARD:

- Do not attempt to open the sterilizer door if a **WATER IN CHAMBER ALARM** condition exists. Call a qualified service technician before attempting to use sterilizer further.
- After manual exhaust, steam may remain inside the chamber. Always wear protective gloves, apron, and a face shield when following emergency procedure to unload sterilizer. Stay as far back from the chamber opening as possible when opening the door.
- Allow sterilizer to cool to room temperature before performing any cleaning or maintenance procedures.
- Failure to shut off the steam supply when cleaning or replacing strainers can result in serious injury.
- Jacket pressure must be 0 psig (0 bar) before beginning work on the steam trap.
- Proper testing of the safety valve requires the valve to be operated under pressure. Exhaust from the safety valve is hot and can cause burns. Proper safety attire (gloves, eye protection, insulated overall) is required. Testing is to be performed by qualified service personnel only.
- When sterilizing liquids, to prevent personal injury or property damage resulting from bursting bottles and hot fluid, you must observe the following procedures:
  - Use **LIQUID** cycle only; no other cycle is safe for processing liquids.
  - Use only vented closures; do not use screw caps or rubber stoppers with crimped seal.
  - Use only Type 1 borosilicate glass bottles; do not use ordinary glass bottles or any container not designed for sterilization.
  - Do not allow hot bottles to be jolted; this can cause hot-bottle explosions. Do not move bottles if any boiling or bubbling is present.



GRAVITY

CANCEL CYCLE START CYCLE



**PURGE PHASE:**

6x (0.41 min)



**CHARGE PHASE:**


45.0 (113.0.F)





Screen #4

**FAST EXHAUST PHASE:**  
 (0.2 inHg) (0.2 inHg) 4.0 (0.2 inHg)




**VACUUM DRY PHASE:**  
 (-0.34 inHg) (-0.34 inHg) 10.0  
 10.0 (-0.34 inHg)



Screen #4

**DRY PHASE:**



Screen #4

AIR BREAK PHASE: 2.0 (-0.17) (-0.07) 2.0 (-0.17)



Screen #4

DEACTIVATE SEAL PHASE: 20



Screen #4

Example of cycle tape.



---

## 6.3 Prevac Cycle (Optional)

---





### WARNING – BURN HAZARD:

- Do not attempt to open the sterilizer door if a WATER IN CHAMBER ALARM condition exists. Call a qualified service technician before attempting to use sterilizer further.
- After manual exhaust, steam may remain inside the chamber. Always wear protective gloves, apron, and a face shield when following emergency procedure to unload sterilizer. Stay as far back from the chamber opening as possible when opening the door.
- Allow sterilizer to cool to room temperature before performing any cleaning or maintenance procedures.
- Failure to shut off the steam supply when cleaning or replacing strainers can result in serious injury.
- Jacket pressure must be 0 psig (0 bar) before beginning work on the steam trap.
- Proper testing of the safety valve requires the valve to be operated under pressure. Exhaust from the safety valve is hot and can cause burns. Proper safety attire (gloves, eye protection, insulated overall) is required. Testing is to be performed by qualified service personnel only.
- When sterilizing liquids, to prevent personal injury or property damage resulting from bursting bottles and hot fluid, you must observe the following procedures:
  - Use LIQUID cycle only; no other cycle is safe for processing liquids.
  - Use only vented closures; do not use screw caps or rubber stoppers with crimped seal.
  - Use only Type 1 borosilicate glass bottles; do not use ordinary glass bottles or any container not designed for sterilization.
  - Do not allow hot bottles to be jolted; this can cause hot-bottle explosions. Do not move bottles if they are boiling or bubbling is present.



WARNING – EXPLOSION HAZARD: This sterilizer is not designed to process flammable compounds.

PREVACCYCLE. PREVAC




ACTIVATE SEAL PHASE: 20





**PURGE PHASE:**

6 \* (0.41)



**PULSE EXHAUST PHASE:**

(0.27)

(113.0°F)

122.5° (253.0°F)

4.0


(0.2)

45.0°




Scree #4

**PULSE CHARGE PHASE:** 26.0 (0.34)




Scree #4

**CHARGE PHASE:** 45.0° (113.0°F)




Scree #4

**STERILIZE PHASE:** 45.0° (113.0°F)



Scree #4

**FAST EXHAUST PHASE:**  
 (0.2). 4.0 (0.2), 4.0



Scree #4

**VACUUM DRY PHASE:**  
 (-0.34) (10.0) 10.0 (-0.34), 10.0




Scree #4

**DRY PHASE:**




Scree #4

AIR BREAK PHASE:  2.0 (-0.17).



**Example of cycle tape.**



## 6.4 Liquid Cycle



### WARNING – BURN HAZARD:

- Do not attempt to open the sterilizer door if a WATER IN CHAMBER ALARM condition exists. Call a qualified service technician before attempting to use sterilizer further.
- After manual exhaust, steam may remain inside the chamber. Always wear protective gloves, apron, and a face shield when following emergency procedure to unload sterilizer. Stay as far back from the chamber opening as possible when opening the door.
- Allow sterilizer to cool to room temperature before performing any cleaning or maintenance procedures.
- Failure to shut off the steam supply when cleaning or replacing strainers can result in serious injury.
- Jacket pressure must be 0 psig (0 bar) before beginning work on the steam trap.
- Proper testing of the safety valve requires the valve to be operated under pressure. Exhaust from the safety valve is hot and can cause burns. Proper safety attire (gloves, eye protection, insulated overall) is required. Testing is to be performed by qualified service personnel only.
- When sterilizing liquids, to prevent personal injury or property damage resulting from bursting bottles and hot fluid, you must observe the following procedures:
  - Use LIQUID cycle only; no other cycle is safe for processing liquids.
  - Use only vented closures; do not use screw caps or rubber stoppers with crimped seal.
  - Use only Type 1 borosilicate glass bottles; do not use ordinary glass bottles or any container not designed for sterilization.
  - Do not allow hot bottles to be jolted; this can cause hot-bottle explosions. Do not move bottles if any boiling or bubbling is present.



WARNING – EXPLOSION HAZARD: This sterilizer is not designed to process flammable compounds.

LIQUID



Screen #96

CANCEL CYCLE  
 START CYCLE  
 CLOSE DOOR(S)



Screen #4

JACKET CHARGE PHASE:  
 (24 .0°F) 120.0°




Screen #4

ACTIVATE SEAL PHASE:  
 20

**PURGE PHASE:**



**STERILIZE PHASE:**








Screen #4d2 T3.125017 98

**SLOW EXHAUST PHASE:**

(0.17), 1.6 / (0.11), 5.0  
 (0.2), 0.6 / (0.04), 5.0 (0.17), 4.2  
 4.0 (0.2) (0.2)



**VAPOR REMOVAL PHASE:**


10





Scree #4

COMPLETE PHASE:



Scree #19



---

## 6.5 Isothermal Cycle (Optional)

---




### WARNING – BURN HAZARD:

- Do not attempt to open the sterilizer door if a WATER IN CHAMBER ALARM condition exists. Call a qualified service technician before attempting to use sterilizer further.
- After manual exhaust, steam may remain inside the chamber. Always wear protective gloves, apron, and a face shield when following emergency procedure to unload sterilizer. Stay as far back from the chamber opening as possible when opening the door.
- Allow sterilizer to cool to room temperature before performing any cleaning or maintenance procedures.
- Failure to shut off the steam supply when cleaning or replacing strainers can result in serious injury.
- Jacket pressure must be 0 psig (0 bar) before beginning work on the steam trap.
- Proper testing of the safety valve requires the valve to be operated under pressure. Exhaust from the safety valve is hot and can cause burns. Proper safety attire (gloves, eye protection, insulated overall) is required. Testing is to be performed by qualified service personnel only.
- When sterilizing liquids, to prevent personal injury or property damage resulting from bursting bottles and hot fluid, you must observe the following procedures:
  - Use LIQUID cycle only; no other cycle is safe for processing liquids.
  - Use only vented closures; do not use screw caps or rubber stoppers with crimped seal.
  - Use only Type 1 borosilicate glass bottles; do not use ordinary glass bottles or any container not designed for sterilization.
  - Do not allow hot bottles to be jolted; this can cause hot-bottle explosions. Do not move bottles if any boiling or bubbling is present.




WARNING – EXPLOSION HAZARD: This sterilizer is not designed to process flammable compounds.

7 110° (172, 230°F).  
**ISOTHERMAL**




Screen #19S

**CANCEL CYCLE**  
**START CYCLE**  
**CLOSE DOOR(S)**



**CHARGE PHASE:**  
45.0° (113.0°F).


**STERILIZE PHASE:**  
45.0° (113.0°F).



Screen #4



Example of cycle tape.





---


## 6.6 Bowie-Dick Warmup Cycle (Only On Prevacuum Sterilizers)

---



### WARNING – BURN HAZARD:

- Do not attempt to open the sterilizer door if a **WATER IN CHAMBER ALARM** condition exists. Call a qualified service technician before attempting to use sterilizer further.
- After manual exhaust, steam may remain inside the chamber. Always wear protective gloves, apron, and a face shield when following emergency procedure to unload sterilizer. Stay as far back from the chamber opening as possible when opening the door.
- Allow sterilizer to cool to room temperature before performing any cleaning or maintenance procedures.
- Failure to shut off the steam supply when cleaning or replacing strainers can result in serious injury.
- Jacket pressure must be 0 psig (0 bar) before beginning work on the steam trap.




B m C |

**Screen #23**




**Screen #4**




**Screen #4**






Screen #4

**PURGE PHASE:**




Screen #4

**CHARGE PHASE:**



Screen #4


**STERILIZE PHASE:**





Scree #4

**FAST EXHAUST PHASE:**  
 (0.2), 4.0 (0.2), 4.0



Scree #4

**VACUUM DRY PHASE:**  
 (-0.34) (10.0), (-0.34), 10.0



Scree #4

**DRY PHASE:**



Scree #4

AIR BREAK PHASE: 2.0 (-0.17), 2.0 (-0.17)




Scree #4

DEACTIVATE SEAL PHASE: 20



Scree #4

COMPLETE PHASE:



Scree #23

---

## 6.7 Bowie-Dick Cycle (Only On Prevacuum Sterilizers)


---



### WARNING – BURN HAZARD:

- Do not attempt to open the sterilizer door if a **WATER IN CHAMBER ALARM** condition exists. Call a qualified service technician before attempting to use sterilizer further.
- After manual exhaust, steam may remain inside the chamber. Always wear protective gloves, apron, and a face shield when following emergency procedure to unload sterilizer. Stay as far back from the chamber opening as possible when opening the door.
- Allow sterilizer to cool to room temperature before performing any cleaning or maintenance procedures.
- Failure to shut off the steam supply when cleaning or replacing strainers can result in serious injury.
- Jacket pressure must be 0 psig (0 bar) before beginning work on the steam trap.
- Proper testing of the safety valve requires the valve to be operated under pressure. Exhaust from the safety valve is hot and can cause burns. Proper safety attire (gloves, eye protection, insulated overall) is required. Testing is to be performed by qualified service personnel only.
- When sterilizing liquids, to prevent personal injury or property damage resulting from bursting bottles and hot fluid, you must observe the following procedures:
  - Use LIQUID cycle only; no other cycle is safe for processing liquids.
  - Use only vented closures; do not use screw caps or rubber stoppers with crimped seal.
  - Use only Type 1 borosilicate glass bottles; do not use ordinary glass bottles or any container not designed for sterilization.
  - Do not allow hot bottles to be jolted; this can cause hot-bottle explosions. Do not move bottles if they are boiling or bubbling intensely.





```

STERIS DARTTEST, 132.0 C, S= 0.03, 1:22:08 PM 4
LOGIN: STERIS 10/1/2004
PULSE #1 PHASE

```

```

STERIS DARTTEST, 132.0 C, S= 0.03 1:22:47 PM 4
LOGIN: STERIS 10/1/2004
PHASE PULSE #1
STATUS EMACULATING 132.0 C
PRESS 381 PSIG
TEMP 90.7 C
PROJECTED CYCLE COMPLETION TIME

```






Screen #4


**PULSE CHARGE PHASE:**

26.0 (0.34)



**CHARGE PHASE:**

45.0° (113.0°F)







---

## 6.8 Leak Test Cycle (Only On Prevacuum Sterilizers)

---

**⚠ WARNING – BURN HAZARD:**

- Do not attempt to open the sterilizer door if a **WATER IN CHAMBER ALARM** condition exists. Call a qualified service technician before attempting to use sterilizer further.
- After manual exhaust, steam may remain inside the chamber. Always wear protective gloves, apron, and a face shield when following emergency procedure to unload sterilizer. Stay as far back from the chamber opening as possible when opening the door.





PULSE CHARGE PHASE:  $\tau_{ch} = \frac{Q_{ch}}{I_{ch}}$



Screen #4

STABILIZE PHASE: 45.0° (113.0°F).



Screen #4


LEAK TEST PHASE: 10



Screen #4


AIR BREAK PHASE: 2.0 (-0.17).





Screen #4

COMPLETE PHASE:




## 6.9 Cycle Abort



Screen #64



#64,  
#2 #1








#10 RIGHT ARROW #36.




Scree #36

LEFT ARROW #10. GRAVITY



Scree #28

PREVAC #2 GRAVITY



Scree #5

LIQUID GRAVITY  
PREVAC #2 #32

#2  
 ARROW ( LEFT ARROW RIGHT  
 #10. EXIT

**CYCLE NAME**  
 CYCLE NAME=GRAVITY1







**OVERTEMP**

**OVERTEMP=6.0°C**



Screen #39

**WHITE**

0.0° (32.0°F)


40.0° (104.0°F)

**ENTER**

**EXIT**

**OVERDRIVE**

**OVERDRIVE=1.5°C**



Screen #40

**WHITE**

0.0° (32.0°F)


40.0° (104.0°F)

**ENTER**

**EXIT**

**UNDERTEMP**

UNDERTEMP=1.0°C



Scree #41

WHITE

0.0° (32.0°F) 40.0° (104.0°F)

ENTER

#41.

EXIT

#2.

**PRINT INTERVAL**

PRINT INT=2 MIN



Scree #42

WHITE

0

ENTER

#42.

EXIT

#2.

DRY TIME

DRY TIME=0:15:00






WHITE F #6 .) EXIT ENTER #2 .)



Screen #69

**INTERLOCK DOORS (double door sterilizers only)**

INTERLOCK = DOOR B  
 1. NO INTERLOCKS -



Screen #61

2. SELECT DOOR TO UNSEAL - sterile side.  
 3. DOOR A UNSEALED - sterile side.  
 4. DOOR B UNSEALED - sterile side.  
 F M #7 .) TOO LONG IN STEP  
 LEFT ARROW RIGHT ARROW  
 EXIT #7 .)

6.13.2 Too Long In Step

**ACTIVATE SEAL**

**ACTIVATE SEAL=1 MIN**

#57.



**WHITE**

#55.

#57.

**EXIT**

**ENTER**


**AIR BREAK**

**AIR BREAK=5 MIN**

#62.

WHITE #5. EXIT ENTER #55.

DEACTIVATE SEAL  
DEACTIVATE SEAL=1 MIN #63.



WHITE #55. #63. EXIT ENTER

EVACUATE  
EVACUATE=30 MIN #60.





Screen #60

WHITE #60. EXIT ENTER #55.

EXHAUST EXHAUST=10 MIN #5




Screen #58

WHITE #5. EXIT ENTER #55.

JACKET CHARGE JACKET CHARGE=60 MIN #77.







Screen #30.1tio is dis 0.12549

2. SPACE 1 (1).  
SPACE 2 (2 12).
3. WHITE #30. ENTER
4. WHITE.
5. ENTER #30.
6. WHITE. M
7. WHITE. M ( )  
#6  
0
- ENTER #30.
- ENTER #12 ( , 4 butto #12.
10. EXIT #7

6.13.7 Machine Number F

M #7,

MACHINE NUMBER



Screen #46

WHITE #46. EXIT ENTER #7.

6.13.8 Default Values F

M #7

DEFAULTVALUES



Screen #9

YES NO #7.

Musical score for 'SILENCE ALARM'. The score consists of multiple staves of music. Tempo markings '250' and '110' are visible. The text 'SILENCE ALARM' is printed below the staves.

**Alarm**

**Descriptio**













**Scree with Alarm**









Musical notation for the Alarm section, showing a few notes on a staff.

Musical notation for the Description section, including a tempo marking '2 (0.14)'.













Musical notation for the Scree with Alarm section, showing a few notes on a staff.










Musical notation for the Scree with Alarm section, showing a few notes on a staff.

Alarm	Descriptio	Scree with Alarm
	 <p style="text-align: right;">Scree #54</p>	
	 <p style="text-align: right;">Scree #49</p>	
	 <p style="text-align: right;">Scree #4</p>	
		

Alarm	Descriptio	Scree with Alarm
	 <p style="text-align: right; border: 1px solid black; border-radius: 10px; padding: 2px;">Scree #4</p>	
	 <p style="text-align: right; border: 1px solid black; border-radius: 10px; padding: 2px;">Scree #4</p>	
	 <p style="text-align: right; border: 1px solid black; border-radius: 10px; padding: 2px;">Scree #4</p>	
	 <p style="text-align: right; border: 1px solid black; border-radius: 10px; padding: 2px;">Scree #4</p>	



Alarm	Descriptio	Scree with Alarm
	 <p style="text-align: right; border: 1px solid black; border-radius: 10px; padding: 2px 10px;">Scree #4</p>	
	 <p style="text-align: right; border: 1px solid black; border-radius: 10px; padding: 2px 10px;">Scree #4</p>	
	 <p style="text-align: right; border: 1px solid black; border-radius: 10px; padding: 2px 10px;">Scree #2</p>	
	 <p style="text-align: right; border: 1px solid black; border-radius: 10px; padding: 2px 10px;">Scree #2</p>	

Alarm	Descriptio	Screes with Alarm
	 <p style="text-align: right; border: 1px solid black; border-radius: 10px; padding: 2px;">Screes #2</p>	
	 <p style="text-align: right; border: 1px solid black; border-radius: 10px; padding: 2px;">Screes #2</p>	
	 <p style="text-align: right;">Screes #2</p>	

## 8.1 Service Mode

#6.1 SERVICE MODE #6




Screen #6



Screen #13

### 8.1.1 Calibration

CALIBRATION



# TEMPERATURE CALIBRATION

TEMPERATURE


#1



NUMERIC

0.1


2 Pa






# TEMPERATURE LOAD VALUES CALIBRATION

## LOAD VALUES




Screen #71



Screen #72


RTD  
 SET  
 EXIT  
 PRESSURE CALIBRATION  
 PRESSURE #1



Screen #74

LOW 0.0  
 #74.  
 PURGE  
 POINT  
 PRESSURE/TEMP  
 SEAL DOOR  
 LOW WHITE  
 0.0  
 30  
 HIGH

(+/- 0.3\*)  
 35.0\*)  
 PRESSURE/TEMP #74.1  
 EVACUATE  
 (+/- 0.3\*) WHITE VAC MUL  
 2 VAC MUL  
 PRESSURE/TEMP  
 AIR BREAK  
 UNSEAL DOOR  
 COMPLETE  
 PRESSURE LOAD VALUES CALIBRATION  
 LOAD VALUES



EXIT


PRINT CALIBRATION VALUES

PRINT ALL CALIBRATION VALUES

DATE	
/30/2004	16 44
IN	
-----	
0.0	12214
24.4	28846
IN	2.040000
-----	
32.0	321
132.0	1644
-----	
32.0	321
132.0	1567
-----	
32.0	330
132.0	1506
-----	



**DEFAULT CALIBRATION VALUES**  
**DEFAULT ALL CALIBRATION**



YES

NO

8.1.2 I/O Test

**I/O TEST**

: : M  
 / /

**Limit switches:**


- 1
- 1 • A
- 2 • B
- 1 • A
- 2 • B

**AC outputs:**

- 1
- 2
- 3
- 4
- 7
- 35 • A
- 36 • B (\*\*)
- 37 • A
- 3 • B



CONTINUE



Screen #26

CONTINUE



Screen #27

CONTINUE



Screen #37

Musical notation with the word "CONTINUE" written in the center.

1. Musical notation  
2. Musical notation  
7. Musical notation  
0.03 \*0.000

## 8.2 Change Printer Paper Roll (Cybertech Printer)

### PAPER LOADING:

- 1.
2. Feed .1
- 3.
- 4.
5. F
6. Feed
7. F .2.
8. #2. #1
9. F .3.
10. Feed
11. Reprint
- 12.



Figure 8-1. Paper Compartment



**Figure 8-2. Paper Feed Positio**




**Figure 8-3. Paper Guide**



---

## 8.4 Change Printer Ribbon (Mylox Printer)

---









3.0  
3.5  
3.6  
3.7

( )  
.  
( )  
/

1 \*  
1 \*  
2 \*

7.0	F	
7.1		6 *
7.2		6 *
7.3	STERIS TEST.	6 *
7.4		6 *

## 9.2 Spare Parts

- 1.
- 2.
- 3.

NOTE: U o STERISa #o,r # a,r o # q # .U o a #o-  
(.r # #2#r b# # a,r,a .

### REPLACEMENT PARTS

Descriptio	Part Number
------------	-------------

#### Cybertech Pri ter

- |      |               |           |
|------|---------------|-----------|
| * 3, | * 3,* (. 5)   | 3 7352-55 |
| * 3, | * .. (*,i# 5) | 3 7352-55 |

#### Mylox Pri ter

- |      |             |            |
|------|-------------|------------|
| * 3, | * 3,* (. 5) | 3 7352-55  |
| * 3, | * .. (. 2)  | 150 2 -440 |

10.1 PLC Specifications

**PLC: MicroLogix® 1200 Controller\***

Catalog number: 1762-24

Age cy Certificatio : 50, 222.2, .142, 1604, 222.2, .213)

Operating temperature: 0, 55, (32 F, 131 F)

Storage temperature: -40, 5, (-40 F, 15 F)

Operating humidity: 5%, 5% ( )

Variables stored in flash memory for permanent storage.

**ANALOG MODULE: 4 CHANNEL INPUT MODULE**

Catalog number: 1762-F4

Resolution : 15, ( )

Repeatability: +/-0.1%

**RTD MODULE: 4 CHANNEL INPUT MODULE**

Catalog number: 1762-4

Resolution : 0.1, (0.1 F)

Repeatability: +/-0.2, (+/- 0.4 F)

Accuracy: +/-0, (+/- 1.62 F)

Maximum drift: +/-0.026, /, (+/- 0.026 F/F)

**RELAY OUTPUT MODULE: 8 AC OUTPUTS**

Catalog number: 1762-

**DISPLAY: PANELVIEW PLUS® 600**

Catalog number: 2711, - 6, 20

**COMMUNICATION MODE:**

Catalog number: 2711, - 22

\* MicroLogix® 1200 Controller, 1762-24, 1762-F4, 1762-4, 2711-6, 20, 2711-22