

TROLLEY MAINTENANCE LOG

Trolley Serial #: _____

Intake HEPA Change: ____/____ Exhaust HEPA Changed: ____/____

Date	Initials	Air Exhaust Pre-Filter	Intake Pre-filter & Hoses Cleaned	Air Intake Pre-Filter	Trolley Certified	HEPA Filters /Decon
------	----------	------------------------	-----------------------------------	-----------------------	-------------------	---------------------

enter date/initials and boxes

		E2wks	E3mo	E12mo	E12mo	PRN

- Trolley Air Exhaust Pre-Filter is changed every 2 weeks
1. With trolley running pre-filter is spritzed w/Sporicidin & bagged-out.
 2. Trolley exhaust area is wiped w/Sporicidin & new filter installed.
 3. Dirty exhaust filter bag is spritzed w/Sporicidin & removed for disposal

- Trolley Rack Hoses are changed every 3 months
1. Hoses are removed & capped with glove – do not create aerosols.
 2. Hoses are spritzed w/Sporicidin & removed to dirty CW for cleaning/decon

Trolley Air Intake Pre-Filter is changed/ sanitized every 12 months
Trolley Certification every 12 months
= 3 yrs/as needed - indicate intake (IT) or
exhaust (EX) HEPA